

兔抗 AIFM2 多克隆抗体

中文名称:	兔抗 AIFM2 多克隆抗体
英文名称:	Anti-AIFM2 rabbit polyclonal antibody
别名:	AMID; PRG3
抗原:	AIFM2
储存:	冷冻 (-20℃)
宿主:	Rabbit
反应种属:	Human, Mouse
相关类别:	一抗
标记物	Unconjugate
克隆类型:	rabbit polyclonal

技术规格

Background:	This gene encodes a flavoprotein oxidoreductase that binds single stranded DNA and is thought to contribute to apoptosis in the presence of bacterial and viral DNA. The expression of this gene is also found to be induced by tumor suppressor protein p53 in colon cancer cells.
Applications:	ELISA, IHC
Name of antibody:	AIFM2
Immunogen:	Synthetic peptide of human AIFM2
Full name:	apoptosis inducing factor, mitochondria associated 2
Synonyms:	AMID; PRG3
SwissProt:	Q9BRQ8

ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer
IHC Recommend dilution:	20-100

