

兔抗 AKIP1 多克隆抗体

中文名称：兔抗 AKIP1 多克隆抗体

英文名称：Anti-AKIP1 rabbit polyclonal antibody

别名：BCA3; C11orf17

相关类别：一抗

抗原：AKIP1

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a nuclear protein that interacts with protein kinase A catalytic subunit, and regulates the effect of the cAMP-dependent protein kinase signaling pathway on the NF-kappa-B activation cascade. Alternatively spliced transcript variants have been described for this gene.
Applications:	ELISA, IHC
Name of antibody:	AKIP1
Immunogen:	Synthetic peptide of human AKIP1
Full name:	A-kinase interacting protein 1
Synonyms:	BCA3; C11orf17
SwissProt:	Q9NQ31

ELISA Recommended dilution:	5000-10000
IHC positive control:	Human ovarian cancer
IHC Recommend dilution:	25-100

