

## 兔抗 ALDH1L1 多克隆抗体

中文名称：兔抗 ALDH1L1 多克隆抗体

英文名称：Anti-ALDH1L1 rabbit polyclonal antibody

别名：aldehyde dehydrogenase 1 family member L1; FDH; FTHFD; 10-fTHF; 10-FTHFDH

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：ALDH1L1

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

|                      |   |
|----------------------|---|
| <b>Background:</b>   | The protein encoded by this gene catalyzes the conversion of 10-formyltetrahydrofolate, nicotinamide adenine dinucleotide phosphate (NADP+), and water to tetrahydrofolate, NADPH, and carbon dioxide. The encoded protein belongs to the aldehyde dehydrogenase family. Loss of function or expression of this gene is associated with decreased apoptosis, increased cell motility, and cancer progression. There is an antisense transcript that overlaps on the opposite strand with this gene locus. Alternative splicing results in multiple transcript variants. |
| <b>Applications:</b> | ELISA, WB   |

|                                    |  |
|------------------------------------|--|
| <b>Name of antibody:</b>           | ALDH1L1  |
| <b>Immunogen:</b>                  | Fusion protein of human ALDH1L1                  |
| <b>Full name:</b>                  | aldehyde dehydrogenase 1 family member L1        |
| <b>Synonyms:</b>                   | FDH; FTHFD; 10-fTHF; 10-FTHFDH                   |
| <b>SwissProt:</b>                  | O75891   |
| <b>ELISA Recommended dilution:</b> | 5000-10000                                       |
| <b>WB Predicted band size:</b>     | 99 kDa   |
| <b>WB Positive control:</b>        | Mouse kidney tissue and Rat liver tissue lysates |
| <b>WB Recommended dilution:</b>    | 200-1000   |

