

兔抗 ANP32A 多克隆抗体

中文名称:	兔抗 ANP32A 多克隆抗体
英文名称:	Anti-ANP32A rabbit polyclonal antibody
别名:	LANP; MAPM; PP32; HPPCn; PHAP1; PHAPI; I1PP2A; C15orf1
抗原:	ANP32A
储存:	冷冻 (-20℃)
宿主:	Rabbit
相关类别:	一抗
反应种属:	Human
标记物:	Unconjugate
克隆类型:	rabbit polyclonal

技术规格

Background:	Implicated in a number of cellular processes, including proliferation, differentiation, caspase-dependent and caspase-independent apoptosis, suppression of transformation (tumor suppressor), inhibition of protein phosphatase 2A, regulation of mRNA trafficking and stability in association with ELAVL1, and inhibition of acetyltransferases as part of the INHAT (inhibitor of histone acetyltransferases) complex. Plays a role in E4F1-mediated transcriptional repression.
Applications:	ELISA, IHC
Name of antibody:	ANP32A
Immunogen:	Full length fusion protein
Full name:	acidic nuclear phosphoprotein 32 family member A
Synonyms:	LANP; MAPM; PP32; HPPCn; PHAP1; PHAPI; I1PP2A; C15o

	rf1
SwissProt:	P39687
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human colorectal cancer and human cervical cancer
IHC Recommend dilution:	25-100



