

兔抗 ANP32C 多克隆抗体

中文名称: 兔抗 ANP32C 多克隆抗体

英文名称: Anti-ANP32C rabbit polyclonal antibody

别名: PP32R1

相关类别: 一抗

储存: 冷冻 (-20℃)

宿主: Rabbit

抗原: ANP32C

反应种属: Human

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

background:	Phosphoprotein 32 (PP32) is a tumor suppressor that can inhibit several types of cancers, including prostate and breast cancers. The protein encoded by this gene is one of at least two proteins that are similar in amino acid sequence to PP32 and are part of the same acidic nuclear phosphoprotein gene family. However, unlike PP32, the encoded protein is tumorigenic. The tumor suppressor function of PP32 has been localized to a 25 amino acid region that is divergent between PP32 and the protein encoded by this gene. This gene does not contain introns.
Applications:	ELISA, IHC
Name of antibody:	ANP32C

Immunogen:	Synthetic peptide of human ANP32C
Full name:	acidic (leucine-rich) nuclear phosphoprotein 32 family, member C
Synonyms :	PP32R1
SwissProt:	O43423
ELISA Recommended dilution:	1000-2000
IHC positive control:	Human thyroid cancer and Human brain
IHC Recommend dilution:	25-100

