

兔抗 APOD 多克隆抗体

中文名称：兔抗 APOD 多克隆抗体

英文名称：Anti-APOD rabbit polyclonal antibody

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：APOD

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a component of high density lipoprotein that has no marked similarity to other apolipoprotein sequences. It has a high degree of homology to plasma retinol-binding protein and other members of the alpha 2 microglobulin protein superfamily of carrier proteins, also known as lipocalins. This glycoprotein is closely associated with the enzyme lecithin:cholesterol acyltransferase - an enzyme involved in lipoprotein metabolism.
Applications:	ELISA, WB, IHC
Name of antibody:	APOD
Immunogen:	Synthetic peptide of human APOD
Full name:	apolipoprotein D
SwissProt:	P05090

ELISA Recommended dilution:	5000-10000
IHC positive control:	Human brain and human lung cancer
IHC Recommend dilution:	20-100
WB Predicted band size:	21 kDa
WB Positive control:	Human liver tissue lysate
WB Recommended dilution:	200-1000





