

## 兔抗 ARFGEF1 多克隆抗体

中文名称：兔抗 ARFGEF1 多克隆抗体

英文名称：Anti-ARFGEF1 rabbit polyclonal antibody

别名：BIG1; P200; ARFGEP1

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：ARFGEF1

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	ADP-ribosylation factors (ARFs) play an important role in intracellular vesicular trafficking. The protein encoded by this gene is involved in the activation of ARFs by accelerating replacement of bound GDP with GTP. It contains a Sec7 domain, which may be responsible for guanine-nucleotide exchange activity and also brefeldin A inhibition.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	ARFGEF1
<b>Immunogen:</b>	Synthetic peptide of human ARFGEF1
<b>Full name:</b>	ADP ribosylation factor guanine nucleotide exchange factor 1

<b>Synonyms:</b>	BIG1; P200; ARFGEP1
<b>SwissProt:</b>	Q9Y6D6
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human thyroid cancer and human tonsil
<b>IHC Recommend dilution:</b>	20-100



