

兔抗 ARSI 多克隆抗体

中文名称：兔抗 ARSI 多克隆抗体

英文名称： Anti-ARSI rabbit polyclonal antibody

别名： SPG66

抗原： ARSI

储存： 冷冻（-20℃） 避光

宿主： Rabbit

反应种属： Human

相关类别： 一抗

标记物： Unconjugate

克隆类型： Unconjugate

技术规格

Background:	ulfatases (EC 3.1.5.6), such as ARSI, hydrolyze sulfate esters from sulfated steroids, carbohydrates, proteoglycans, and glycolipids. They are involved in hormone biosynthesis, modulation of cell signaling, and degradation of macromolecules (Sardiello et al., 2005 [PubMed 16174644])
Applications:	WB, IF
Name of antibody:	ARSI
Immunogen:	Synthesized peptide derived from internal of human ARSI.
Full name:	arylsulfatase family, member I
Synonyms :	SPG66

SwissProt:	Q5FYB1
WB Predicted band size:	64 kDa
WB Positive control:	COS7 cells lysate
WB Recommended dilution:	500-3000
IF positive control:	MCF-7 cells
IF Recommend dilution:	100-500



