

兔抗 ASTE1 多克隆抗体

中文名称: 兔抗 ASTE1 多克隆抗体

英文名称: Anti-ASTE1 rabbit polyclonal antibody

别名: asteroid homolog 1; HT001

储存: 冷冻 (-20℃)

抗原: ASTE1

宿主: Rabbit

相关类别: 一抗

反应种属: Human

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	Possible role in EGF receptor signaling.
Applications:	ELISA, WB, IHC
Name of antibody:	ASTE1
Immunogen:	Fusion protein of human ASTE1
Full name:	asteroid homolog 1
Synonyms:	HT001
SwissProt:	Q2TB18
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer
IHC Recommend dilution:	50-300
WB Predicted band size:	77 kDa
WB Positive control:	A172 cell lysate

WB Recommended dilution:

200-1000



