

兔抗 ATAD1 多克隆抗体

中文名称：兔抗 ATAD1 多克隆抗体

英文名称：Anti-ATAD1 rabbit polyclonal antibody

别名：AFDC1; FNP001; THORASE

相关类别：一抗

抗原：ATAD1

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	ATPase that plays a critical role in regulating the surface expression of AMPA receptors (AMPA), thereby regulating synaptic plasticity and learning and memory. Required for NMDA-stimulated AMPAR internalization and inhibition of GRIA1 and GRIA2 recycling back to the plasma membrane; these activities are ATPase-dependent (By similarity).
Applications:	ELISA, WB, IHC
Name of antibody:	ATAD1
Immunogen:	Fusion protein of human ATAD1
Full name:	ATPase family, AAA domain containing 1
Synonyms:	AFDC1; FNP001; THORASE

SwissProt:	Q8NBU5
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer and human esophagus cancer
IHC Recommend dilution:	50-300
WB Predicted band size:	41 kDa
WB Positive control:	Rat brain tissue, Mouse brain tissue and Human cerebrum tissue lysates
WB Recommended dilution:	500-2000





