

兔抗 ATAD3A 多克隆抗体

中文名称：兔抗 ATAD3A 多克隆抗体

英文名称：Anti-ATAD3A rabbit polyclonal antibody

相关类别：一抗

储 存： 冷冻（-20℃）

宿 主： Rabbit

抗 原： ATAD3A

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

| | |
|------------------------------------|--|
| Background: | ATAD3A and ATAD3B (MIM 612317) are mitochondrial membrane proteins that contribute to the stabilization of large mitochondrial DNA (mtDNA)-protein complexes called nucleoids (He et al., 2007 [PubMed 17210950]). |
| Applications: | ELISA, WB, IHC |
| Name of antibody: | ATAD3A |
| Immunogen: | Fusion protein of human ATAD3A |
| Full name: | ATPase family, AAA domain containing 3A |
| SwissProt: | Q9NVI7 |
| ELISA Recommended dilution: | 5000-10000 |
| IHC positive control: | Human lung cancer and human prostate cancer |
| IHC Recommend dilution: | 25-100 |
| WB Predicted band size: | 71 kDa |

| | |
|---------------------------------|---------------------------------|
| WB Positive control: | Raji,Hela and A431 cell lysates |
| WB Recommended dilution: | 500-2000 |





