

兔抗 **BATF** 多克隆抗体

中文名称: 兔抗 BATF 多克隆抗体

英文名称: Anti-BATF rabbit polyclonal antibody

别名: SFA2; B-ATF; BATF1; SFA-2

储存: 冷冻 (-20℃)

抗原: BATF

宿主: Rabbit

相关类别: 一抗

反应种属: Human, Mouse, Rat

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is a nuclear basic leucine zipper protein that belongs to the AP-1/ATF superfamily of transcription factors. The leucine zipper of this protein mediates dimerization with members of the Jun family of proteins. This protein is thought to be a negative regulator of AP-1/ATF transcriptional events.
Applications:	ELISA, IHC
Name of antibody:	BATF
Immunogen:	Full length fusion protein
Full name:	basic leucine zipper ATF-like transcription factor
Synonyms:	SFA2; B-ATF; BATF1; SFA-2
SwissProt:	Q16520

ELISA Recommended dilution:	5000-10000
IHC positive control:	Human thyroid cancer and human ovarian cancer
IHC Recommend dilution:	30-150



