

## 兔抗 BICD2 多克隆抗体

中文名称：兔抗 BICD2 多克隆抗体

英文名称：Anti-BICD2 rabbit polyclonal antibody

别名：SMALED2; bA526D8.1

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：BICD2

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene is one of two human homologs of Drosophila bicaudal-D and a member of the Bicoid family. It has been implicated in dynein-mediated, minus end-directed motility along microtubules. It has also been reported to be a phosphorylation target of NIMA related kinase 8. Two alternative splice variants have been described.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	BICD2
<b>Immunogen:</b>	Synthetic peptide of human BICD2
<b>Full name:</b>	BICD cargo adaptor 2
<b>Synonyms:</b>	SMALED2; SMALED2A; SMALED2B; bA526D8.1
<b>SwissProt:</b>	Q8TD16

<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human liver cancer
<b>IHC Recommend dilution:</b>	20-100

