

兔抗 BLM 多克隆抗体

中文名称: 兔抗 BLM 多克隆抗体

英文名称: Anti-BLM rabbit polyclonal antibody

别名: BS; RECQ2; RECQL2; RECQL3

相关类别: 一抗

储存: 冷冻 (-20℃)

宿主: Rabbit

抗原: BLM

反应种属: Human

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	The Bloom syndrome gene product is related to the RecQ subset of DExH box-containing DNA helicases and has both DNA-stimulated ATPase and ATP-dependent DNA helicase activities. Mutations causing Bloom syndrome delete or alter helicase motifs and may disable the 3'-5' helicase activity. The normal protein may act to suppress inappropriate recombination.
Applications:	ELISA, IHC
Name of antibody:	BLM
Immunogen:	Synthetic peptide of human BLM
Full name:	Bloom syndrome RecQ like helicase
Synonyms:	BS; RECQ2; RECQL2; RECQL3

SwissProt:	P54132
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human colorectal cancer and human tonsil
IHC Recommend dilution:	25-100

