

## 兔抗 BLVRA 多克隆抗体

中文名称：兔抗 BLVRA 多克隆抗体

英文名称：Anti-BLVRA rabbit polyclonal antibody

别名：biliverdin reductase A; BVR; BLVR; BVRA

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：BLVRA

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	The protein encoded by this gene belongs to the biliverdin reductase family, members of which catalyze the conversion of biliverdin to bilirubin in the presence of NADPH or NADH. Mutations in this gene are associated with hyperbiliverdinemia. Alternatively spliced transcript variants have been found for this gene.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	BLVRA
<b>Immunogen:</b>	Full length fusion protein
<b>Full name:</b>	biliverdin reductase A
<b>Synonyms:</b>	BVR; BLVR; BVRA
<b>SwissProt:</b>	P53004

<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human gastric cancer and Human colorectal cancer
<b>IHC Recommend dilution:</b>	20-100



