

兔抗 BMP2 多克隆抗体

中文名称：兔抗 BMP2 多克隆抗体

英文名称：Anti-BMP2 rabbit polyclonal antibody

别名：BDA2; BMP2A

储存：冷冻 (-20℃)

抗原：BMP2

宿主：Rabbit

相关类别：一抗

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a member of the transforming growth factor-beta (TGFB) superfamily. The encoded preproprotein is proteolytically processed to generate each subunit of the disulfide-linked homodimer, which induces bone and cartilage formation.
Applications:	ELISA, WB, IHC
Name of antibody:	BMP2
Immunogen:	Synthetic peptide of human BMP2
Full name:	bone morphogenetic protein 2
Synonyms:	BDA2; BMP2A
SwissProt:	P12643
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer
IHC Recommend dilution:	40-200

WB Predicted band size:	45 kDa
WB Positive control:	Human fetal liver tissue lysate
WB Recommended dilution:	500-2000

