

兔抗 BRD8 多克隆抗体

中文名称: 兔抗 BRD8 多克隆抗体

英文名称: Anti-BRD8 rabbit polyclonal antibody

别名: SMAP; p120; SMAP2

储存: 冷冻 (-20℃)

抗原: BRD8

宿主: Rabbit

相关类别: 一抗

反应种属: Human, Mouse

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	The protein encoded by this gene interacts with thyroid hormone receptor in a ligand-dependent manner and enhances thyroid hormone-dependent activation from thyroid response elements. This protein contains a bromodomain and is thought to be a nuclear receptor coactivator. Multiple alternatively spliced transcript variants that encode distinct isoforms have been identified.
Applications:	ELISA, IHC
Name of antibody:	BRD8
Immunogen:	Synthetic peptide of human BRD8
Full name:	bromodomain containing 8

Synonyms:	SMAP; p120; SMAP2
SwissProt:	Q9H0E9
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human esophagus cancer and human ovarian cancer
IHC Recommend dilution:	25-100

