

兔抗 BRF1 多克隆抗体

中文名称: 兔抗 BRF1 多克隆抗体

英文名称: Anti-BRF1 rabbit polyclonal antibody

别名: BRF; CFDS; hBRF; BRF-1; GTF3B; TAF3C; TAF3B2; TF3B90; TAFIII90; TFIIB90; HEL-S-76p

抗原: BRF1

储存: 冷冻 (-20℃)

宿主: Rabbit

反应种属: Human

相关类别: 一抗

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	This gene encodes one of the three subunits of the RNA polymerase III transcription factor complex. This complex plays a central role in transcription initiation by RNA polymerase III on genes encoding tRNA, 5S rRNA, and other small structural RNAs. The gene product belongs to the TF2B family. Several alternatively spliced variants encoding different isoforms, that function at different promoters transcribed by RNA polymerase III, have been identified.
Applications:	ELISA, IHC
Name of antibody:	BRF1

Immunogen:	Full length fusion protein
Full name:	BRF1, RNA polymerase III transcription initiation factor 90 kDa subunit
Synonyms:	BRF; CFDS; hBRF; BRF-1; GTF3B; TAF3C; TAF3B2; TF3B90; TAFIII90; TFIIIB90; HEL-S-76p
SwissProt:	Q92994
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human ovarian cancer and human thyroid cancer
IHC Recommend dilution:	25-100



