

兔抗 CALB2 多克隆抗体

中文名称: 兔抗 CALB2 多克隆抗体

英文名称: Anti-CALB2 rabbit polyclonal antibody

别名: CR; CAL2; CAB29

抗原: CALB2

储存: 冷冻 (-20℃)

宿主: Rabbit

相关类别: 一抗

反应种属: Human, Mouse, Rat

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	This gene encodes an intracellular calcium-binding protein belonging to the troponin C superfamily. Members of this protein family have six EF-hand domains which bind calcium. This protein plays a role in diverse cellular functions, including message targeting and intracellular calcium buffering. It also functions as a modulator of neuronal excitability, and is a diagnostic marker for some human diseases, including Hirschsprung disease and some cancers. Alternative splicing results in multiple transcript variants.
Applications:	ELISA, WB, IHC
Name of antibody:	CALB2

Immunogen:	Synthetic peptide of human CALB2
Full name:	calbindin 2
Synonyms :	CR; CAL2; CAB29
SwissProt:	P22676
ELISA Recommended dilution:	1000-2000
IHC positive control:	Human colorectal cancer
IHC Recommend dilution:	10-50
WB Predicted band size:	32 kDa
WB Positive control:	293T cell
WB Recommended dilution:	200-1000



