

兔抗 CASS4 多克隆抗体

中文名称：兔抗 CASS4 多克隆抗体

英文名称：Anti-CASS4 rabbit polyclonal antibody

别名：Cas scaffolding protein family member 4; CAS4; HEFL; HEPL; C20orf32

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：CASS4

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Docking protein that plays a role in tyrosine kinase-based signaling related to cell adhesion and cell spreading. Regulates PTK2/FAK1 activity, focal adhesion integrity, and cell spreading.
Applications:	ELISA, IHC
Name of antibody:	CASS4
Immunogen:	Synthetic peptide of human CASS4
Full name:	Cas scaffolding protein family member 4
Synonyms:	CAS4; HEFL; HEPL; C20orf32
SwissProt:	Q9NQ75
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human tonsil

IHC Recommend dilution:

20-100

