

兔抗 CBX7 多克隆抗体

中文名称: 兔抗 CBX7 多克隆抗体

英文名称: Anti-CBX7 rabbit polyclonal antibody

储 存: 冷冻 (-20℃)

抗 原: CBX7

宿 主: Rabbit

反应种属: Human, Mouse, Rat

相关类别: 一抗

标 记 物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	CBX7 (Chromobox protein homolog 7) is a 251 amino acid nuclear protein that contains one N-terminal chromo domain and one C-terminal Pc box. Highly expressed in kidney, brain, heart and skeletal muscle, with weaker expression in peripheral blood leukocytes, CBX7 functions as a component of the chromatin-associated polycomb complex (PcG) and is involved in maintaining the transcriptionally repressed state of target genes.
Applications:	ELISA, WB
Name of antibody:	CBX7
Immunogen:	Synthetic peptide of human CBX7
Full name:	chromobox homolog 7
SwissProt:	O95931
ELISA Recommended dilution:	2000-5000

WB Predicted band size:	28 kDa
WB Positive control:	Hela, Jurkat, 293T and lovo cell
WB Recommended dilution:	500-2000

