

## 兔抗 CCDC85C 多克隆抗体

中文名称: 兔抗 CCDC85C 多克隆抗体

英文名称: Anti-CCDC85C rabbit polyclonal antibody

抗 原: CCDC85C

储 存: 冷冻 (-20℃)

宿 主: Rabbit

反应种属: Human, Mouse

相关类别: 一抗

标 记 物: Unconjugate

克隆类型: rabbit polyclonal

### 技术规格

<b>Background:</b>	CCDC85C may play an important role in cortical development, especially in the maintenance of radial glia.
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	CCDC85C
<b>Immunogen:</b>	Synthetic peptide of human CCDC85C
<b>Full name:</b>	coiled-coil domain containing 85C
<b>SwissProt:</b>	A6NKD9
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human liver cancer
<b>IHC Recommend dilution:</b>	20-100
<b>WB Predicted band size:</b>	45 kDa
<b>WB Positive control:</b>	Human breast cancer tissue lysate
<b>WB Recommended dilution:</b>	200-1000



