

兔抗 CCNB1IP1 多克隆抗体

中文名称：兔抗 CCNB1IP1 多克隆抗体

英文名称：Anti-CCNB1IP1 rabbit polyclonal antibody

别名：cyclin B1 interacting protein 1; HEI10; C14orf18

相关类别：一抗

储存：冷冻（-20℃）

抗原：CCNB1IP1

宿主：Rabbit

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	HEI10 is a member of the E3 ubiquitin ligase family and functions in progression of the cell cycle through G(2)/M.[supplied by OMIM, Apr 2004]
Applications:	ELISA, IHC
Name of antibody:	CCNB1IP1
Immunogen:	Full length fusion protein
Full name:	cyclin B1 interacting protein 1
Synonyms:	HEI10; C14orf18
SwissProt:	Q9NPC3
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human lung cancer and human colorectal cancer
IHC Recommend dilution:	30-150



