

## 兔抗 CDC14A 多克隆抗体

中文名称：兔抗 CDC14A 多克隆抗体

英文名称：Anti-CDC14A rabbit polyclonal antibody

别名：cdc14; hCDC14

相关类别：一抗

储存：冷冻（-20℃）避光

宿主：Rabbit

抗原：CDC14A

反应种属：Human

标记物：Unconjugate

克隆类型：Unconjugate

### 技术规格

<b>Background:</b>	The protein encoded by this gene is a member of the dual specificity protein tyrosine phosphatase family. It is highly similar to <i>Saccharomyces cerevisiae</i> Cdc14, a protein tyrosine phosphatase involved in the exit of cell mitosis and initiation of DNA replication, suggesting a role in cell cycle control. This protein has been shown to interact with, and dephosphorylate tumor suppressor protein p53, and is thought to regulate the function of p53. Alternative splicing of this gene results in several transcript variants encoding distinct isoforms.
<b>Applications:</b>	WB
<b>Name of antibody:</b>	CDC14A

<b>Immunogen:</b>	Synthesized peptide derived from Internal of human C C14A.
<b>Full name:</b>	cell division cycle 14A
<b>Synonyms :</b>	cdc14; hCDC14
<b>SwissProt:</b>	Q9UNH5
<b>WB Predicted band size:</b>	67 kDa
<b>WB Positive control:</b>	Hela cells lysate
<b>WB Recommended dilution:</b>	500-3000

