

兔抗 CDK20 多克隆抗体

中文名称：兔抗 CDK20 多克隆抗体

英文名称：Anti-CDK20 rabbit polyclonal antibody

别名：cyclin dependent kinase 20; P42; CCRK; CDCH; PNQALRE

相关类别：一抗

抗原：CDK20

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The protein encoded by this gene contains a kinase domain most closely related to the cyclin-dependent protein kinases. The encoded kinase may activate cyclin-dependent kinase 2 and is involved in cell growth. Alternatively spliced transcript variants encoding distinct isoforms have been reported.
Applications:	ELISA, IHC
Name of antibody:	CDK20
Immunogen:	Synthetic peptide of human CDK20
Full name:	cyclin dependent kinase 20
Synonyms:	P42; CCRK; CDCH; PNQALRE
SwissProt:	Q8IZL9

ELISA Recommended dilution:	1000-2000
IHC positive control:	Human liver cancer and human lung cancer
IHC Recommend dilution:	20-100



