

兔抗 CDV3 多克隆抗体

中文名称：兔抗 CDV3 多克隆抗体

英文名称：Anti-CDV3 rabbit polyclonal antibody

别名：CDV3 homolog; H41

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：CDV3

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The CDV3 protein is a 258 amino acid protein that localizes in the cytoplasm. It acts as a repressor of N eu and subsequently may be involved in the progression of some cancers, including breast cancer, gliomas and prostate tumors.
Applications:	ELISA, IHC
Name of antibody:	CDV3
Immunogen:	Synthetic peptide of human CDV3
Full name:	CDV3 homolog
Synonyms:	H41
SwissProt:	Q9UKY7
ELISA Recommended dilution:	5000-10000

IHC positive control:	Human liver cancer
IHC Recommend dilution:	10-50

