

## 兔抗 CELF5 多克隆抗体

中文名称：兔抗 CELF5 多克隆抗体

英文名称：Anti-CELF5 rabbit polyclonal antibody

别名：CELF-5; BRUNOL5; BRUNOL-5

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：CELF5

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a member of the the CELF/BRUNOL protein family, which contain two N-terminal RNA recognition motif (RRM) domains, one C-terminal RRM domain, and a divergent segment of 160-230 aa between the second and third RRM domains. Members of this protein family regulate pre-mRNA alternative splicing and may also be involved in mRNA editing and translation. Alternatively spliced transcript variants have been found for this gene.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	CELF5
<b>Immunogen:</b>	Synthetic peptide of human CELF5

<b>Full name:</b>	CUGBP, Elav-like family member 5
<b>Synonyms:</b>	CELF-5; BRUNOL5; BRUNOL-5
<b>SwissProt:</b>	Q8N6W0
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human tonsil and human thyroid cancer
<b>IHC Recommend dilution:</b>	40-200



