

兔抗 CHD4 多克隆抗体

中文名称：兔抗 CHD4 多克隆抗体

英文名称： Anti-CHD4 rabbit polyclonal antibody

别名： chromodomain helicase DNA binding protein 4

相关类别： 一抗

储存： 冷冻（-20℃） 避光

抗原： CHD4

宿主： Rabbit

反应种属： Human Mouse

标记物： Unconjugate

克隆类型： Unconjugate

技术规格

Background:	The product of this gene belongs to the SNF2/RAD54 helicase family. It represents the main component of the nucleosome remodeling and deacetylase complex and plays an important role in epigenetic transcriptional repression. Patients with dermatomyositis develop antibodies against this protein. Somatic mutations in this gene are associated with serous endometrial tumors. Alternative splicing results in multiple transcript variants encoding different isoforms.
Applications:	IHC, IF
Name of antibody:	CHD4
Immunogen:	Synthesized peptide derived from internal of human CHD4.

Full name:	chromodomain helicase DNA binding protein 4
Synonyms :	chromodomain helicase DNA binding protein 4
SwissProt:	Q14839
IHC positive control:	Human liver carcinoma tissue
IHC Recommend dilution:	50-100
IF positive control:	HepG2 cells
IF Recommend dilution:	100-500



