

兔抗 CHFR 多克隆抗体

中文名称： 兔抗 CHFR 多克隆抗体

英文名称： Anti-CHFR rabbit polyclonal antibody

别名： RNF116; RNF196

相关类别： 一抗

储存： 冷冻（-20℃） 避光

宿主： Rabbit

抗原： CHFR

反应种属： Human Mouse Rat

标记物： Unconjugate

技术规格

Background:	This gene encodes an E3 ubiquitin-protein ligase required for the maintenance of the antephasal checkpoint that regulates cell cycle entry into mitosis and, therefore, may play a key role in cell cycle progression and tumorigenesis. The encoded protein has an N-terminal forkhead-associated domain, a central RING-finger domain, and a cysteine-rich C-terminal region. Alternatively spliced transcript variants that encode different protein isoforms have been described.
Applications:	WB
Name of antibody:	CHFR
Immunogen:	Synthesized peptide derived from Internal of human CHFR.
Full name:	checkpoint with forkhead and ring finger domains, E3 ubiquitin protein ligase

Synonyms :	RNF116; RNF196
SwissProt:	Q96EP1
WB Predicted band size:	73 kDa
WB Positive control:	Jurkat cells lysate
WB Recommended dilution:	500-3000

