

兔抗 CKS2 多克隆抗体

- 中文名称：兔抗 CKS2 多克隆抗体
- 英文名称：Anti-CKS2 rabbit polyclonal antibody
- 别名：CDC28 protein kinase regulatory subunit 2; CKSHS2
- 相关类别：一抗
- 抗原：CKS2
- 储存：冷冻（-20℃）
- 宿主：Rabbit
- 反应种属：Human, Mouse
- 标记物：Unconjugate
- 克隆类型：rabbit polyclonal

技术规格

Background:	CKS2 protein binds to the catalytic subunit of the cyclin dependent kinases and is essential for their biological function. The CKS2 mRNA is found to be expressed in different patterns through the cell cycle in HeLa cells, which reflects specialized role for the encoded protein.
Applications:	ELISA, IHC
Name of antibody:	CKS2
Immunogen:	Synthetic peptide of human CKS2
Full name:	CDC28 protein kinase regulatory subunit 2
Synonyms:	CKSHS2

SwissProt:	P33552
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human thyroid cancer
IHC Recommend dilution:	25-100

