

兔抗 CLIC6 多克隆抗体

中文名称：兔抗 CLIC6 多克隆抗体

英文名称： Anti-CLIC6 rabbit polyclonal antibody

别名： chloride intracellular channel 6; CLIC1L

相关类别： 一抗

抗原： CLIC6

储存： 冷冻（-20℃）

宿主： Rabbit

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes a member of the chloride intracellular channel family of proteins. The gene is part of a large triplicated region found on chromosomes 1, 6, and 21. Alternative splicing results in multiple transcript variants encoding different isoforms.
Applications:	ELISA, IHC
Name of antibody:	CLIC6
Immunogen:	Synthetic peptide of human CLIC6
Full name:	chloride intracellular channel 6
Synonyms:	CLIC1L
SwissProt:	Q96NY7

ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer
IHC Recommend dilution:	30-150



