

## 兔抗 CLTCL1 多克隆抗体

中文名称：兔抗 CLTCL1 多克隆抗体

英文名称： Anti-CLTCL1 rabbit polyclonal antibody

别名： clathrin heavy chain like 1; CLTD; CHC22; CLH22; CLTCL

相关类别： 一抗

抗原： CLTCL1

储存： 冷冻（-20℃）

宿主： Rabbit

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene is a member of the clathrin heavy chain family and encodes a major protein of the polyhedral coat of coated pits and vesicles. Chromosomal aberrations involving this gene are associated with meningioma, Di George syndrome, and velo-cardio-facial syndrome. Multiple transcript variants encoding different isoforms have been found for this gene.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	CLTCL1
<b>Immunogen:</b>	Synthetic peptide of human CLTCL1
<b>Full name:</b>	clathrin heavy chain like 1

<b>Synonyms:</b>	CLTD; CHC22; CLH22; CLTCL
<b>SwissProt:</b>	P53675
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human liver cancer and human ovarian cancer
<b>IHC Recommend dilution:</b>	40-200



