

兔抗 COL8A2 多克隆抗体

中文名称: 兔抗 COL8A2 多克隆抗体

英文名称: Anti-COL8A2 rabbit polyclonal antibody

别名: FECD; PPCD; FECD1; PPCD2

抗原: COL8A2

储存: 冷冻 (-20℃)

宿主: Rabbit

反应种属: Human, Mouse

相关类别: 一抗

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	This gene encodes the alpha 2 chain of type VIII collagen. This protein is a major component of the basement membrane of the corneal endothelium and forms homo- or heterotrimers with alpha 1 (VIII) type collagen. Defects in this gene are associated with Fuchs endothelial corneal dystrophy and posterior polymorphous corneal dystrophy type 2. Alternative splicing results in multiple transcript variants.
Applications:	ELISA, IHC
Name of antibody:	COL8A2
Immunogen:	Synthetic peptide of human COL8A2
Full name:	collagen, type VIII, alpha 2

Synonyms :	FECD; PPCD; FECD1; PPCD2
SwissProt:	P25067
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human liver cancer and Human esophagus cancer
IHC Recommend dilution:	25-100

