

兔抗 COLEC12 多克隆抗体

中文名称: 兔抗 COLEC12 多克隆抗体

英文名称: Anti-COLEC12 rabbit polyclonal antibody

别名: CLP1; NSR2; SRCL; SCARA4

抗原: COLEC12

储存: 冷冻 (-20℃)

宿主: Rabbit

相关类别: 一抗

反应种属: Human, Mouse, Rat

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	This gene encodes a member of the C-lectin family, proteins that possess collagen-like sequences and carbohydrate recognition domains. This protein is a scavenger receptor, a cell surface glycoprotein that displays several functions associated with host defense. It can bind to carbohydrate antigens on microorganisms, facilitating their recognition and removal. It also mediates the recognition, internalization, and degradation of oxidatively modified low density lipoprotein by vascular endothelial cells.
Applications:	ELISA, IHC
Name of antibody:	COLEC12
Immunogen:	Synthetic peptide of human COLEC12

Full name:	collectin sub-family member 12
Synonyms :	CLP1; NSR2; SRCL; SCARA4
SwissProt:	Q5KU26
ELISA Recommended dilution:	2000-10000
IHC positive control:	Human brain and Human esophagus cancer
IHC Recommend dilution:	30-150

