

## 兔抗 CORO2B 多克隆抗体

中文名称: 兔抗 CORO2B 多克隆抗体

英文名称: Anti-CORO2B rabbit polyclonal antibody

别名: coronin 2B; CLIPINC

相关类别: 一抗

储存: 冷冻 (-20℃)

宿主: Rabbit

抗原: CORO2B

反应种属: Human, Mouse

标记物: Unconjugate

克隆类型: rabbit polyclonal

### 技术规格

<b>Background:</b>	May play a role in the reorganization of neuronal actin structure.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	CORO2B
<b>Immunogen:</b>	Synthetic peptide of human CORO2B
<b>Full name:</b>	coronin 2B
<b>Synonyms:</b>	CLIPINC
<b>SwissProt:</b>	Q9UQ03
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human gastric cancer
<b>IHC Recommend dilution:</b>	25-100

