

兔抗 CPLX1 多克隆抗体

中文名称: 兔抗 CPLX1 多克隆抗体

英文名称: Anti-CPLX1 rabbit polyclonal antibody

别名: complexin 1; CPX1; CPX-I

抗原: CPLX1

储存: 冷冻 (-20℃)

宿主: Rabbit

反应种属: Human, Mouse, Rat

相关类别: 一抗

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	Proteins encoded by the complexin/synaphin gene family are cytosolic proteins that function in synaptic vesicle exocytosis. These proteins bind syntaxin, part of the SNAP receptor. The protein product of this gene binds to the SNAP receptor complex and disrupts it, allowing transmitter release.
Applications:	ELISA, IHC
Name of antibody:	CPLX1
Immunogen:	Fusion protein of human CPLX1
Full name:	complexin 1
Synonyms:	CPX1; CPX-I
SwissProt:	O14810

ELISA Recommended dilution:	5000-10000
IHC positive control:	Human cervical cancer and Human thyroid cancer
IHC Recommend dilution:	25-100



