

## 兔抗 CPSF2 多克隆抗体

- 中文名称: 兔抗 CPSF2 多克隆抗体
- 英文名称: Anti-CPSF2 rabbit polyclonal antibody
- 别名: cleavage and polyadenylation specific factor 2 CPSF100
- 相关类别: 一抗
- 储存: 冷冻 (-20℃)
- 宿主: Rabbit
- 抗原: CPSF2
- 反应种属: Human, Mouse
- 标记物: Unconjugate
- 克隆类型: rabbit polyclonal

### 技术规格

<b>Background:</b>	Component of the cleavage and polyadenylation specificity factor (CPSF) complex that play a key role in pre-mRNA 3'-end formation, recognizing the AAUAAA signal sequence and interacting with poly(A) polymerase and other factors to bring about cleavage and poly(A) addition. Involved in the histone 3' end pre-mRNA processing.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	CPSF2
<b>Immunogen:</b>	Synthetic peptide of human CPSF2
<b>Full name:</b>	cleavage and polyadenylation specific factor 2
<b>Synonyms:</b>	CPSF100

<b>SwissProt:</b>	Q9P2I0
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human esophagus cancer
<b>IHC Recommend dilution:</b>	25-100

