

## 兔抗 CR1L 多克隆抗体

中文名称：兔抗 CR1L 多克隆抗体

英文名称：Anti-CR1L rabbit polyclonal antibody

别名：complement C3b/C4b receptor 1 like

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：CR1L

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	The C3b/C4b receptor (CR1; 120620) is important in immune complex processing. On phagocytic cells it promotes the adherence and sometimes the endocytosis of C3b- and C4b-coated particles. CR1 is composed largely of a tandemly repeated motif of approximately 60 amino acids, the short consensus repeat (SCR).
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	CR1L
<b>Immunogen:</b>	Synthetic peptide of human CR1L
<b>Full name:</b>	complement C3b/C4b receptor 1 like
<b>SwissProt:</b>	Q2VPA4
<b>ELISA Recommended dilution:</b>	5000-10000

<b>IHC positive control:</b>	Human liver cancer and Human gastric cancer
<b>IHC Recommend dilution:</b>	20-100



