

## 兔抗 CREB3L2 多克隆抗体

- 中文名称: 兔抗 CREB3L2 多克隆抗体
- 英文名称: Anti-CREB3L2 rabbit polyclonal antibody
- 别名: cAMP responsive element binding protein 3 like 2 BBF2H7
- 储存: 冷冻 (-20℃)
- 抗原: CREB3L2
- 宿主: Rabbit
- 反应种属: Human
- 标记物: Unconjugate
- 相关类别: 一抗
- 克隆类型: rabbit polyclonal

### 技术规格

<b>Background:</b>	Full length fusion protein
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	CREB3L2
<b>Immunogen:</b>	Full length fusion protein
<b>Full name:</b>	cAMP responsive element binding protein 3 like 2
<b>Synonyms:</b>	BBF2H7
<b>SwissProt:</b>	Q70SY1
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human liver cancer;
<b>IHC Recommend dilution:</b>	25-100
<b>WB Predicted band size:</b>	57 kDa
<b>WB Positive control:</b>	HEPG2 and Hela cell lysates

WB Recommended dilution: 200-1000



