

## 兔抗 CSNK1A11L 多克隆抗体

中文名称：兔抗 CSNK1A1/1L 多克隆抗体

英文名称： Anti-CSNK1A1/1L rabbit polyclonal antibody

别名： CK1; CK1a; CK1 $\alpha$ ; HLCDGP1; PRO2975; HEL-S-77p/CK1

相关类别： 一抗

储存： 冷冻（-20℃） 避光

抗原： CSNK1A1/1L

宿主： Rabbit

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

|                              |  |
|------------------------------|--|
| <b>Background:</b>           | Casein kinases are operationally defined by their preferential utilization of acidic proteins such as caseins as substrates. It can phosphorylate a large number of proteins. Participates in Wnt signaling By similarity. |
| <b>Applications:</b>         | IHC, IF  |
| <b>Name of antibody:</b>     | CSNK1A1/1L   |
| <b>Immunogen:</b>            | Synthesized peptide derived from internal of human CKI- $\alpha$ 1/1-like  |
| <b>Full name:</b>            | casein kinase 1, alpha 1/1-like  |
| <b>Synonyms :</b>            | CK1; CK1a; CK1 $\alpha$ ; HLCDGP1; PRO2975; HEL-S-77p/CK1  |
| <b>SwissProt:</b>            | Q8N752/P48729  |
| <b>IHC positive control:</b> | Human testis tissue  |

|                                |            |
|--------------------------------|------------|
| <b>IHC Recommend dilution:</b> | 50-100     |
| <b>IF positive control:</b>    | A549 cells |
| <b>IF Recommend dilution:</b>  | 100-500    |



