

兔抗 CS 多克隆抗体

中文名称: 兔抗 CS 多克隆抗体

英文名称: Anti-CS rabbit polyclonal antibody

相关类别: 一抗

储 存: 冷冻 (-20℃)

宿 主: Rabbit

抗 原: CS

反应种属: Human, Mouse, Rat

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is a Krebs tricarboxylic acid cycle enzyme that catalyzes the synthesis of citrate from oxaloacetate and acetyl coenzyme A. The enzyme is found in nearly all cells capable of oxidative metabolism. This protein is nuclear encoded and transported into the mitochondrial matrix, where the mature form is found.
Applications:	ELISA, WB, IHC
Name of antibody:	CS
Immunogen:	Fusion protein of human CS
Full name:	citrate synthase
SwissProt:	O75390
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human ovarian cancer

IHC Recommend dilution:	25-100
WB Predicted band size:	52 kDa
WB Positive control:	Jurkat and 293T cell lysates
WB Recommended dilution:	500-2000



