

兔抗 CUBN 多克隆抗体

中文名称：兔抗 CUBN 多克隆抗体

英文名称：Anti-CUBN rabbit polyclonal antibody

别名：IFCR; MGA1; gp280

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：CUBN

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：Unconjugate

技术规格

Background:	Cubilin (CUBN) acts as a receptor for intrinsic factor-vitamin B12 complexes. The role of receptor is supported by the presence of 27 CUB domains. Cubilin is located within the epithelium of intestine and kidney. Mutations in CUBN may play a role in autosomal recessive megaloblastic anemia.
Applications:	ELISA, WB, IHC
Name of antibody:	CUBN
Immunogen:	Synthetic peptide of human CUBN
Full name:	cubilin (intrinsic factor-cobalamin receptor)
Synonyms :	IFCR; MGA1; gp280
SwissProt:	O60494

ELISA Recommended dilution:	5000-10000
IHC positive control:	Human tonsil
IHC Recommend dilution:	30-150
WB Predicted band size:	399 kDa
WB Positive control:	Human kidney tissue
WB Recommended dilution:	200-1000



