

## 兔抗 CWF19L1 多克隆抗体

中文名称：兔抗 CWF19L1 多克隆抗体

英文名称：Anti-CWF19L1 rabbit polyclonal antibody

别名：CWF19-like 1, cell cycle control (S. pombe) C19L1; hDrn1; SCAR17

抗原：CWF19L1

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human, Mouse

标记物：Unconjugate

相关类别：一抗

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a member of the CWF19 protein family. Mutations in this gene have been associated with autosomal recessive spinocerebellar ataxia-17 and mild cognitive disability. Alternative splicing results in multiple transcript variants.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	CWF19L1
<b>Immunogen:</b>	Fusion protein of human CWF19L1
<b>Full name:</b>	CWF19 like cell cycle control factor 1
<b>Synonyms:</b>	C19L1; hDrn1; SCAR17
<b>SwissProt:</b>	Q69YN2
<b>ELISA Recommended dilution:</b>	5000-10000

<b>IHC positive control:</b>	Human prostate cancer
<b>IHC Recommend dilution:</b>	20-100

