

## 兔抗 CX3CR1 多克隆抗体

中文名称: 兔抗 CX3CR1 多克隆抗体

英文名称: Anti-CX3CR1 rabbit polyclonal antibody

别名: V28; CCRL1; GPR13; CMKDR1; GPRV28; CMKBRL1

相关类别: 一抗

抗原: CX3CR1

储存: 冷冻 (-20℃)

宿主: Rabbit

反应种属: Human

标记物: Unconjugate

克隆类型: Unconjugate

### 技术规格

<b>Background:</b>	Fractalkine is a transmembrane protein and chemokine involved in the adhesion and migration of leukocytes. The protein encoded by this gene is a receptor for fractalkine. The encoded protein also is a coreceptor for HIV-1, and some variations in this gene lead to increased susceptibility to HIV-1 infection and rapid progression to AIDS. Four transcript variants encoding two different isoforms have been found for this gene.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	CX3CR1
<b>Immunogen:</b>	Synthetic peptide of human CX3CR1
<b>Full name:</b>	chemokine (C-X3-C motif) receptor 1

<b>Synonyms :</b>	V28; CCRL1; GPR13; CMKDR1; GPRV28; CMKBRL1
<b>SwissProt:</b>	P49238
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human lung cancer
<b>IHC Recommend dilution:</b>	25-100

