

兔抗 CYP20A1 多克隆抗体

中文名称：兔抗 CYP20A1 多克隆抗体

英文名称：Anti-CYP20A1 rabbit polyclonal antibody

别名：CYP-M

相关类别：一抗

储存：冷冻（-20℃）避光

宿主：Rabbit

抗原：CYP20A1

反应种属：Human

标记物：Unconjugate

克隆类型：Unconjugate

技术规格

Background:	This gene encodes a member of the cytochrome P450 superfamily of enzymes. The cytochrome P450 proteins are monooxygenases that catalyze many reactions involved in drug metabolism and synthesis of cholesterol, steroids and other lipids. This protein lacks one amino acid of the conserved heme binding site. It also lacks the conserved I-helix motif A GX(D,E)T, suggesting that its substrate may carry its own oxygen.
Applications:	WB
Name of antibody:	CYP20A1
Immunogen:	Synthesized peptide derived from internal of human Cytochrome P450 20A1.

Full name:	cytochrome P450, family 20, subfamily A, polypeptide 1
Synonyms :	CYP-M
SwissProt:	Q6UW02
WB Predicted band size:	52 kDa
WB Positive control:	Jurkat cells lysate
WB Recommended dilution:	500-3000

