

兔抗 CYP4X1 多克隆抗体

中文名称：兔抗 CYP4X1 多克隆抗体

英文名称：Anti-CYP4X1 rabbit polyclonal antibody

别名：CYP4X1

相关类别：一抗

储存：冷冻（-20℃）避光

宿主：Rabbit

抗原：CYP4X1

反应种属：Human

标记物：Unconjugate

克隆类型：Unconjugate

技术规格

Background:	This gene encodes a member of the cytochrome P450 superfamily of enzymes. The cytochrome P450 proteins are monooxygenases which catalyze many reactions involved in drug metabolism and synthesis of cholesterol, steroids and other lipids. The expression pattern of a similar rat protein suggests that this protein may be involved in neurovascular function in the brain.
Applications:	WB
Name of antibody:	CYP4X1
Immunogen:	Synthesized peptide derived from internal of human Cytochrome P450 4X1.
Full name:	cytochrome P450, family 4, subfamily X, polypeptide 1

Synonyms :	CYPIVX1
SwissProt:	Q8N118
WB Predicted band size:	59 kDa
WB Positive control:	293 cells and Jurkat cells lysates
WB Recommended dilution:	500-3000

