

## 兔抗 DDIT3 多克隆抗体

中文名称: 兔抗 DDIT3 多克隆抗体

英文名称: Anti-DDIT3 rabbit polyclonal antibody

别名: CHOP; CEBPZ; CHOP10; CHOP-10; GADD153

相关类别: 一抗

储存: 冷冻 (-20℃)

宿主: Rabbit

抗原: DDIT3

反应种属: Human, Mouse, Rat

标记物: Unconjugate

克隆类型: Unconjugate

### 技术规格

**Background:**

This gene encodes a member of the CCAAT/enhancer-binding protein (C/EBP) family of transcription factors. The protein functions as a dominant-negative inhibitor by forming heterodimers with other C/EBP members, such as C/EBP and LAP (liver activator protein), and preventing their DNA binding activity. The protein is implicated in adipogenesis and erythropoiesis, is activated by endoplasmic reticulum stress, and promotes apoptosis. Fusion of this gene and FUS on chromosome 16 or EWSR1 on chromosome 22 induced by translocation generates chimeric proteins in myxoid liposarcomas or Ewing sarcoma. Multiple alternatively spliced transcript variants encoding t

	wo isoforms with different length have been identified.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	DDIT3
<b>Immunogen:</b>	Synthetic peptide of human DDIT3
<b>Full name:</b>	DNA-damage-inducible transcript 3
<b>Synonyms :</b>	CHOP; CEBPZ; CHOP10; CHOP-10; GADD153
<b>SwissProt:</b>	P35638
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human tonsil
<b>IHC Recommend dilution:</b>	25-100

