

## 兔抗 DEF8 多克隆抗体

中文名称: 兔抗 DEF8 多克隆抗体

英文名称: Anti-DEF8 rabbit polyclonal antibody

别名: differentially expressed in FDCP 8 homolog

相关类别: 一抗

抗原: DEF8

储存: 冷冻 (-20℃)

宿主: Rabbit

反应种属: Human

标记物: Unconjugate

克隆类型: rabbit polyclonal

### 技术规格

<b>Background:</b>	The function of this protein remains unknown.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	DEF8
<b>Immunogen:</b>	Full length fusion protein
<b>Full name:</b>	differentially expressed in FDCP 8 homolog
<b>SwissProt:</b>	Q6ZN54
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human thyroid cancer
<b>IHC Recommend dilution:</b>	30-150

