

兔抗 DEFA5 多克隆抗体

中文名称：兔抗 DEFA5 多克隆抗体

英文名称：Anti-DEFA5 rabbit polyclonal antibody

别名：DEF5; HD-5

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：DEFA5

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Defensins are a family of microbicidal and cytotoxic peptides thought to be involved in host defense. They are abundant in the granules of neutrophils and also found in the epithelia of mucosal surfaces such as those of the intestine, respiratory tract, urinary tract, and vagina. Members of the defensin family are highly similar in protein sequence and distinguished by a conserved cysteine motif. Several of the alpha defensin genes appear to be clustered on chromosome 8. The protein encoded by this gene, defensin, alpha 5, is highly expressed in the secretory granules of Paneth cells of the ileum.
Applications:	ELISA, IHC

Name of antibody:	DEFA5
Immunogen:	Synthetic peptide of human DEFA5
Full name:	defensin, alpha 5, Paneth cell-specific
Synonyms:	DEF5; HD-5
SwissProt:	Q01523
IHC positive control:	Human gastric cancer and human brain
IHC Recommend dilution:	15-50

