

兔抗 DEFB112 多克隆抗体

中文名称: 兔抗 DEFB112 多克隆抗体

英文名称: Anti-DEFB112 rabbit polyclonal antibody

别名: DEFB-12

抗原: DEFB112

储存: 冷冻 (-20℃)

宿主: Rabbit

相关类别: 一抗

反应种属: Human

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	Defensins form a family of antimicrobial and cytotoxic peptides made by neutrophils. Defensins are short, processed peptide molecules that are classified by structure into three groups: alpha-defensins, beta-defensins and theta-defensins. All beta-defensin genes are densely clustered in four to five syntenic chromosomal regions.
Applications:	ELISA, IHC
Name of antibody:	DEFB112
Immunogen:	Synthetic peptide of human DEFB112
Full name:	defensin, beta 112

Synonyms :	DEFB-12
SwissProt:	Q30KQ8
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human ovarian cancer
IHC Recommend dilution:	20-100

